

Recommended Magnification Size

Application	Magnification Range
General Surgery	2.5x-3.0x
Spine Surgery	3.0x-3.5x
Neurosurgery	3.5x-4.0x
Plastic and reconstructive Surgery	3.5x-5.0x
Cardiac Surgery	3.0x-3.5x
Ophthalmic surgery	3.0x-3.5x

2.5x magnification loupes are used by general surgeons to obtain the correct posture and avoid eye strain. 6x magnification is used for delicate micro surgery procedures. Loupe selection is more of a personal preference and comfort. The table above gives you the guiding recommendations for different surgical areas.

STAC Loupes are backed by excellent customer service of ACTS Medical and a world wide network of qualified local representatives.

Discuss your needs and preferences with our professional representatives. They will help you to select the right product for your needs.



All STAC Loupes come with

Side Shields, Carrying Bag, Spanner with Fork, Cleaning Cloth and Screwdriver as Standard Accessories.

ACTS Medical

STAC Microsurgery Loupes



Surgery Made Easier

ACTS Medical



A Product of ACTS Medical, Canada

STAC Loupes

Your Friend in Deed



As per recent study 93% surgeons suffer eye strain and 67% suffer lower back

Many micro surgeons suffer a shortened professional career due to occupational hazards such as Eye Strain and lower back pain.. STAC Loupes are specially designed to improve visual acuity and help prevent eye strain and lower back and neck pain. STAC loupes optics are aligned using state of the art laser technology and go through rigorous testing to give you better image visualization, enhanced depth of field and wide range of working distance. Whether you are working on Craniotomy procedure, Lymph node extirpation procedure, Coronary artery by-

pass grafting procedure, Bisegmental cervical spinal canal stenosis procedure or any other microsurgery procedure, STAC Loupes are your friend in deed.



Lymph node extirpation, sub-mandibular left, for chronically enlarged lymph nodes

STAC Standard Features

- ⇒ Laser aligned optics for crisp focus
- ⇒ Titanium frame for lighter weight and sturdiness
- ⇒ Choice of 4 working distances for perfect working position
- ⇒ Choice of 6 models, magnification power from 2.5x to 6x to meet all your surgical needs
- ⇒ Few simple steps to change your prescription glasses .
- ⇒ One step adjustment to intrapapillary convergence and declination angle (up to 45°) for perfect optical alignment and body posture..

STAC Advantage

- ⇒ Your patients will benefit by a higher standard of quality care
- ⇒ You will benefit by enhanced vision, body-friendly ergonomics, less fatigue, and an extended career.



Actual Size



3.5x Magnification

Analytical & Medical Technologies



SINCE 2006

Analyticals Medical Technologies Pvt.Ltd.
Shree Ambika Heritage, 1st Floor, Plot No.1, Sector 1,
Near Kharghar Railway Station, Navi Mumbai 410210
www.amtmedicals.com
info@amtmedicals.com

+91 8976074090